# ImmunoCAP™ Allergen

## Fluoroenzymimmunoassay

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## USO PRETENDIDO

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#### REFERÊNCIAS

- Wide L, Bennich H, Johansson SGO. Diagnosis of allergy by an in vitro test for allergen antibodies. *Lancet* 1967;2:1105-7.
- Johansson SGO, Yman L. In vitro assays for immunoglobulin E. *Clin Rev Allergy* 1988;6(2):93-139.
- Akdis M et al.: Mechanisms of allergen specific immunotherapy. *J Allergy Clin Immunol* 2007;119:780-789.
- Francis JN, Grass pollen immunotherapy: IL-10 induction and suppression of late responses precedes IgG4 inhibitory antibody activity. *J Allergy Clin Immunol* 2008;121:1120-5.
- Lima MT, Wilson D, Pitkin L et al. Grass pollen sublingual immunotherapy for seasonal rhinoconjunctivitis: a randomized controlled trial. *Clin Exp Allergy* 2002;32:507-14.
- Ewbank PA, Murray J, Sanders K, Curran-Everett D, Dreskin S, Nelson HS. A double-blind, placebo-controlled immunotherapy dose-response study with standardized cat extract. *J Allergy Clin Immunol* 2003; 111:155-61.
- Rossi RE et al. Evaluation of serum IgG4 antibodies specific to grass pollen allergen components in the follow up of allergic patients undergoing subcutaneous and sublingual immunotherapy. *Vaccine* 2007; 25:957-964.
- C. Bussmann et al. Clinical improvement and immunological changes in atopic dermatitis patients undergoing subcutaneous immunotherapy with a house dust mite allergoid: a pilot study *Clinical and Experimental Allergy*, 2007; 37:1277-1285.
- Stapel SO, Asero R, Ballmer-Weber BK, Knol EF, Strobel S., Vieths S., Kleine- Tebbe J.. Testing for IgG4 against foods is not recommended as a diagnostic tool: EAACI Task Force Report. *Allergy* 2008; 63: 793-796.

- Platts-Mills T, Vaughan J, Squillace S, Woodfolk J, Sporik R. Sensitisation, asthma and a modified Th2 response in children exposed to cat allergen: a populationbased cross-sectional study. *Lancet* 2001; 357:752-6.
- Lemon-Mule H et al. Immunologic changes in children with egg allergy ingesting extensively heated egg *J Allergy Clin Immunol* 2008;122(5):977- 983.e1. Epub 2008 Oct 11.
- Akdis M. Healthy immune response to allergens: T regulatory cells and more. *Curr Opin Immunol* 2006;18:738-744.
- Moverare R, Vestermen E, Metsot, Sorva R, Elfman L, Haahtela T: Pollenspecific rush immunotherapy: clinical efficacy and effects on antibody concentrations. *Ann Allergy Asthma Immunol* 2001;86(3):337-42.
- Ljorring C, Frew AJ, Durham S, Jacobsen L, Johansen N, Nedergaard- Larsen Jet al.: Induction of a Competitive Non-IgE Antibody Response by Specific Immunotherapy is Related to Clinical Efficacy. *J Allergy Clin Immunol* 2006;117(2):264.
- Buchanan A, Green TD, Jones SM, Scurlock AM, Christie L, Althage KA, Steele PH, Pons L, Helm RM, Lee LA, and Wesley Burks A: Egg oral immunotherapy in nonanaphylactic children with egg allergy *J Allergy Clin Immunol* 2007;119:199-205.
- Matsui EC, Diette GB, Kropw EJM, Aalberse RC, Smith AL, Curtin-Brosnan J, Eggleston PA: Mouse allergen-specific immunoglobulin G and immunoglobulin G4 and allergic symptoms in immunoglobulin E-sensitized laboratory animal workers. *Clin Exp Allergy* 2005;35:1347-1353.
- Beyer K, Teuber SS: Food allergy diagnostics: scientific and unproven procedures. *Curr Opin Allergy Clin Immunol*. 2005 Jun;5(3):261-6.
- van Hoeyveld E, Dupont L, Bossuyt X: Quantification of IgG antibodies to Aspergillus fumigatus and Pigeon antigens by ImmunoCAP technology: An alternative to the precipitation technique? *Clin Chem* 2006;52:1785-1793.
- Barton RC, Hobson RP, Denton M, Peckham D, Brownlee K, Conway S, Kerr MA: Serologic diagnosis of allergic bronchopulmonary aspergillosis in patients with cystic fibrosis through the detection of immunoglobulin G to Aspergillus fumigatus. *Diagn Microbiol Infect Dis*. 2008 Nov;62(3):287-91.
- Rydjord B, Hetland G, Wiker HG: Immunoglobulin G antibodies against environmental moulds in a Norwegian healthy population shows a bimodal distribution for Aspergillus versicolor. *Scand J Immunol*. 2005 Sep;62(3):281-8.
- CSLI. Analytical Performance Characteristics and Clinical Utility of Immunological Assays for Human Immunoglobulin E (IgE) Antibodies and Defined Allergen Specificities; Approved Guideline – Second Edition. I/LA20-A2 Vol. 29 No. 9.
- Tate J, Ward G: Interference in Immunoassay. *Clin Biochem Rev* 2004;25:105-120.
- Guttormsen A.B, Johansson S.G.O, Öman H., Wilhelmssen V., Nopp A. No consumption of IgE antibody in serum during allergic drug anaphylaxis. *Allergy* 2007; 62: 1326-1330.
- Goldberg A., Confino-Cohen R. Timing of venom skin tests and IgE determinations after insect sting anaphylaxis. *J Allergy Clin Immunol* 1997; 100:182-184.
- Wide L, Juhlin L. Detection of penicillin allergy of the immediate type by radio-immunoassay of reagins (IgE) to penicilloyl conjugates. *Clin Allergy* 1971; 1: 171-7.
- Kraft D, Wide L. Clinical patterns and results of radioallergosorbent test (RAST) and skin tests in penicillin allergy. *Br J Dermatol* 1976;94:593-601.

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© Phadia AB, Uppsala, Sweden.
Phadia AB,