## Calprotectin 2

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## REFERENCIAS

1. Gaya DR, Mackenzie JF (2002). Fecal calprotectin: a bright future for assessing disease activity in Crohn's disease. *Q J Med* 95: 557-558
2. Roseth AG et al (2004). Normalization of fecal calprotectin: a predictor of mucosal healing in patients with inflammatory bowel disease. *Scand J Gastroenterol* 39: 1017-1020
3. Vermeire S et al (2006). Laboratory markers in IBD: useful, magic or unnecessary toys? *Gut* 55: 426-431
4. Masoodi I et al (2011). Biomarkers in the management of ulcerative colitis: a brief review. *Ger Med Sci* 9: Doc03. doi: 10.3205/000126
5. Tibble J et al (2000). A simple method for assessing intestinal inflammation in Crohn's disease. *Gut* 47: 506-513
6. Sutherland AD et al (2008). Review of fecal biomarkers in inflammatory bowel disease. *Dis Colon Rectum* 51:1283-1291
7. Poullis A et al (2004). Bowel inflammation as measured by fecal calprotectin: A link between lifestyle factors and colorectal cancer risk. *Cancer Epidemiol Biomarkers Prev* 13: 279-284
8. Roca M et al (2017). Fecal calprotectin and eosinophil-derived neurotoxin in healthy children between 0-12 years. *J Pediatr Gastroenterol Nutr.* 2017 Oct;65(4):394-398

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© Phadia GmbH, Freiburg, Germany

## Revision History

Version	Countries	Change
21	all	Phadia 250 and Phadia 2500/5000, chapter "Summary and explanation of the test": Insertion of the sentence "Fecal calprotectin levels get affected by nonsteroidal anti-inflammatory drug (NSAID) intake, bleeding more than 100 ml and by malignancy".
		Phadia 250 and Phadia 2500/5000, chapter "Specimen collection, handling and preparation", at the end of both parts A and B: Incorporation of the stool extract stability.
		Phadia 250 and Phadia 2500/5000, chapter "Limitations": Deletion of the sentence "Blood in the stool may affect EliA Calprotectin results".
		Phadia 250 and Phadia 2500/5000, chapter "References": Insertion of a new literature article (reference number 4) and renumbering of the following references.
22	all	Harmonization of texts and translations between DfUs. For more details, please contact your local representative.
		Chapter "Reagents/Material": for Stop Solution (Art. No. 34-2337-11), changed storage temperatures to 2–32°C.
		Chapter "Expected values": rephrasing of first paragraph to "The normal population tests positive for calprotectin in a certain percentage, which is strongly dependent on age. While it has been reported that 11% of children are positive, the corresponding number for adults between 50 and 70 years of age is 25%."
22	UK	Chapter "Specimen collection, handling and preparation", part A: change of storage conditions for the extract from 4°C to 2–8°C.
22	CZ	Chapter "Reagents/Material": EliA Calprotectin Well, addition of "Ready for use" and "dry" to storage conditions.
		Chapter "Reagents/Material": Calprotectin Positive Control, "Lidský calprotectinu v PBS obsahujícím BSA..." changed to "Lidský calprotectin v PBS obsahujícím BSA...".
		Chapter "Reagents/Material": EliA Sample Diluent (Art. No. 83-1023-01), corrected bottle volume from 64 ml to 48 ml.
		Chapter "Reagents/Material": EliA Stool Extraction Kit 2 (Art. No. 14-5651-01), corrected buffer volume to 1300 µl.
		Chapter "Reagents/Material": Stop Solution (Art. No. 10-9442-01), correction of number of bottles from 4 to 6.
		Chapter "Warnings and precautions": insertion of "Some of the reagents are manufactured from human blood components".
		Chapter "Specimen collection, handling and preparation": fill volume for extraction tubes corrected to 1300 µl.
		Chapter "Specimen collection, handling and preparation": extract storage conditions: Polish text replaced with Czech.
		Chapter "Specimen collection, handling and preparation": addition of transfer volume 1-2 ml to section B-7.

22	CZ	Chapter "Procedure" for Phadia 250: correction of Well name to Calprotectin and addition of extra sentence "In case of manual handling make sure to enter the characters below the barcode."
		Chapter "Procedure": updated texts on Handling of EliA Calprotectin 2 Well.
		Chapter "Procedure", on-board stability of reagents: addition of missing storage conditions for EliA Wells for Phadia 250: "...or 24 hours at room temperature."
		Chapter "Procedure", on-board stability of reagents for Phadia 250: addition of EliA Calibrator Strips, EliA Curve Control Strips and their storage conditions.
		Chapter "Procedure", on-board stability of reagents, EliA Sample Diluent and Stop Solution for Phadia 250: addition of "Re-cap bottles every night".
		Chapter "Procedure", on-board stability of reagents, EliA Conjugate for Phadia 250: corrected information.
		Chapter "Procedure", on-board stability of reagents, Development Solution for Phadia 250: corrected on-board stability conditions.
		Chapter "Procedure", on-board stability of reagents, Washing Solution for Phadia 250: additional information added.
		Chapter "Procedure", on-board stability of reagents: "Washing Solution Concentrate/Additive" deleted from Phadia 250.
		Chapter "Procedure", Procedural comments for Phadia 250: Addition of "From one sample diluted by the instrument (1:200), up to 11 determinations can be made."
		Chapter "Procedure", Procedural comments for Phadia 2500/5000: updated information about availability of first results.
		Chapter "Calibration and reference material": corrected text on International standards for Calprotectin.
22	DE/AT/CH	Harmonization of the three German language DfUs (German, Austrian and Swiss) into one DfU for all three languages.
22	DE	Chapter "Reagents/Material": addition of "Positive" to Control name: "Alle Komponenten außer der EliA Calprotectin Positive Control ..."
22	ES	Chapter "Procedure": for Phadia 2500/5000, table note under Sample volume per determination table: corrected text "sin volumen de cebado. Para el volumen de cebado ver el Manual de Usuario Phadia 2500/5000."
22	FR	Chapter "Reagents/Material": Correction of buffer from "Calprotectine humaine dans une solution de tampon Tris contenant ..." to "Calprotectine humaine dans une solution de tampon phosphate contenant ..."
22	IT	Chapter "Specimen collection, handling and preparation": upper storage limit for extract corrected to ≤ -20 °C.
		Chapter "Performance characteristics", Specificity: addition of "human" to "la determinazione della calprotectina umana".
22	NO	Chapter "Reagents/Material": EliA Calprotectin Well, addition of "dry" to storage conditions.

22	PL	Chapter "Reagents/Material": Development Solution (Art. No. 10-9440-01) correction of number of determinations to ">170"
		Chapter "Procedure", On-board stability for Phadia 250: Norwegian text replaced with Polish.
22	PT	Chapter "Reagents/Material": for Phadia 250, EliA Calprotectin Positive Control 100 (Art. No. 83-1066-01) changed to EliA Calprotectin Positive Control 250 (Art. No. 83-1083-01).
		Chapter "Reagents/Material": for Phadia 250, note under Positive Control, EliA Calprotectin Positive Control 100 changed to EliA Calprotectin Positive Control 250.
		Chapter "Reagents/Material": for Phadia 250, EliA Calprotectin Negative Control 100 (Art. No. 83-1067-01) changed to EliA Calprotectin Negative Control 250 (Art. No. 83-1085-01).
		Chapter "Reagents/Material": for Phadia 250, note under Negative Control, EliA Calprotectin Positive Control 100 changed to EliA Calprotectin Negative Control 250.
22	SE	Chapter "Principles of the procedure": text correction "När ej bundna antikroppar tvättats bort tillsätts enzymmärkta antikroppar mot humana Calprotectin-antikroppar (EliA Calprotectin 2-konjugat) för att bilda ett calprotectinkonjugat-komplex." changed to "När ej bundna komponenter tvättats bort tillsätts enzymmärkta antikroppar mot humant calprotectin (EliA Calprotectin 2-konjugat) för att bilda ett calprotectin-konjugat-komplex."
		Chapter "Reagents/Material": Stop Solution (Art. no. 34-2337-01), "50 rör (18 wells vardera)" corrected to "6 flaskor (2800 ml vardera); reagenser för >13100 bestämningar".
		Chapter "Reagents/Material": Dilution Wells (Art. no. 12-4005-69), "60 plattor per förpackning; tillräckligt för 60 x 96 prover" corrected to "50 rör (18 brunnar vardera)".
		Chapter "Reagents/Material": addition of "wedge-shaped" to EliA Calprotectin 2 Conjugate 50 and 200.
22	NO, PL, SE	Chapter "Reagents/Material": addition of "wedge-shaped" to EliA Calprotectin 2 Conjugate 50 and 200.
22	IT, PL	Chapter "Reagents/Material": Stop Solution (Art. No. 10-9442-01), correction of number of determinations to ">560"
22	CZ, PL	Chapter "Reagents/Material": EliA Calprotectin 2 Extraction Buffer (Art. No. 83-1147-01), corrected volume from 115 ml to 117 ml.
		Chapter "Reagents/Material": EliA Calprotectin Calibrator Strips and EliA Calprotectin 2 Curve Control Strips, corrected unit from ng/l to ng/ml.
		Chapters "Calculation and presentation of results" and "Performance characteristics": corrected units from ng/l to ng/ml.
22	CZ, DE, ES, FR, IT, PT, SE	Chapter "Reagents/Material": for Phadia 250 General Reagents, addition of Development Solution (Art. No. 10-9441-01) and Stop Solution (Art. No. 10-9479-01)
22	CZ, NO, PL	Chapter "Reagents/Material": Washing Solution for Phadia 250, deletion of additional information regarding shelf life and materials not supplied.

22	CZ, DK, ES, NO, PT, SE	Chapter "Calibration and reference material": "There are no international standards for calprotectin. Results are given in arbitrary mg/kg", deletion of "arbitrary".
22	CZ, DE, DK, ES, FR, IT, NO, PL, PT, SE	Chapter "Calculation and presentation of results": "Phadia 250 measures specific calprotectin concentrations in ng/ml", deletion of "specific".
23	all, except US	Harmonization of texts and translations between DfUs/instruments.
		General text improvements throughout the DfU, e.g. changing "sample(s)" to "extract(s)".
		Chapter "Reagents/Material": rephrasing the composition of EliA Calprotectin Conjugates.
		Phadia 200 and Phadia 250, chapter "Reagents/Material", Dilution Plates: Remove information on plastic material.
		Phadia 200, chapter "Reagents/Material", Dilution Plates: Correction of package size from 50 to 100.
		Chapter "Warnings and Precautions": insertion of "Do not use if desiccant bag is missing or foilbag is damaged."
		Chapter "Warnings and Precautions": insertion of "Wear gloves while handling samples and reagents provided."
		Phadia 200, Phadia 250 and Phadia 2500/5000, chapter "Instrument": Clarification that the Phadia instrument processes all steps of the test.
		Chapter "Specimen Collection, Handling and Extract Preparation", update of extraction method A.
		Phadia 250 and Phadia 2500/5000, chapter "Procedure": Update and new presentation of On-board stability section.
		Phadia 200, chapter "Procedure", On-board stability, Development Solution and Stop Solution: Addition of "Recap bottles every night".
		Chapter "Expected Values": Updated text.
		Chapter "Performance Characteristics", Specificity: Updated text.
Phadia 200, chapter "Performance Characteristics", Precision: corrected explanation of study design.		
Deletion of reference number 7, renumbering of references after that. New literature article reference 8 added.		
Updated symbol table. Link to full symbol glossary inserted.		
23	CZ	Correction of Czech translations.
23	BG	Chapter "Reagents/Material": Correction of translations.
23	NL	Chapters "Reagents/Material", "Calculation and Interpretation of Results" and "Performance Characteristics": Updated unit from ng/l to ng/ml.
23	GR	Chapter "Specimen Collection, Handling and Extract Preparation", Extract Dilution: Correction of translations.