



ABON™ One-Step Ovulation Tests

Empowering Pregnancy Planning

For The Right Time To Conceive!

- Accurate
- Fast
- Easy to Use



ABON™ One-Step Ovulation Tests

Empowering Pregnancy Planning

Accurate

- ≥99% compared to other lab tests
- Predicts ovulation within 24-36 hours
- Available in options of 40 mIU/mL or 30 mIU/mL sensitivity

Fast

- 3-5 minutes



Easy to Use

- Add Specimen → Wait 3-5 minutes → Read Results

Test Procedure


Strip Test

- Immerse the strip into urine
- Read results at 5 minutes





Device Test

- Add 3 drops of urine
- Read results at 5 minutes



Midstream Test

- Pass the urine on the absorbent tip or dip it into urine
- Read results at 3-5 minutes



Ordering Information

Product Name	Specimen	Format	Catalog No.	Sensitivity	Tests per Kit
LH Ovulation	Urine	Strip	FLH-101 † †	40 mIU/mL	50
		Device	FLH-102 † †	40 mIU/mL	40
			FLH-102H †	40 mIU/mL	5
		Midstream	FLH-103 † †	40 mIU/mL	5
			FLH-U103 †	30 mIU/mL	

† US FDA Cleared † CE Marked CE



Inverness Medical Innovations Hong Kong Limited
 #15th Floor, Building B, Modern International Plaza,
 159 Tianmu Shan Road, Hangzhou 310007, P.R.China
 Please visit our website for details: www.imihongkong.com
 E-mail: IMIHK.info@invmed.com

TECNIGEN S.A.
 Av. Zañartu 1060
 Ñuñoa, Santiago
 Chile
www.tecnigen.cl
info@tecnigen.cl

